Mui Scientific

Division of
H&A Mui
Enterprises Inc.

145 Traders Blvd. E., Unit #34, Mississauga, ON, Canada, L4Z 3L3
Tel: (905) 890-5525 Toll Free: 1-800-303-6611 Fax: (905) 890-3523
Website: www.muiscientific.com Email: mail@muiscientific.com

REF CR1B SERIES: R MODEL: 1B

EN: REUSABLE PVC MANOMETRIC CATHETER
ANORECTAL EXPLUSION CATHETER (AIR ONLY)
INSTRUCTIONS FOR USE

INTENDED PURPOSE

The Anorectal Exulsion Catheter is a specific Reusable PVC Manometric Catheter model which isolates the anorectal balloon function of an anorectal manometric catheter. The Anorectal Exulsion Catheter is used to evaluate the muscular contractions of the rectum.

INDICATIONS FOR USE

The use of the Anorectal Exulsion Catheter is indicated when there is a need or suspicion of an anorectal disorder.

CONTRAINDICATIONS

The presence of any stenosis, structural deformity or disease within the rectum, that results in an unacceptably increased risk of distention, perforation, or any other adverse consequence arising from use of this device.

WARNINGS

Read through entire IFU to reduce any possible risks from misuse. This catheter is supplied non-sterile. Inspect catheter and packaging for any evidence of mechanical damage or imperfections. Do not use if damaged. This catheter must not be altered or modified in any way. Mui Scientific is not liable for personal injury and damage to property if original Mui Scientific parts are not being used.

Self-impaction of catheter balloon is performed at your own risk. Proper aseptic technique and universal barrier precautions (UBP) must apply. Post-procedure manometric catheters are contaminated with body fluid. Immediately clean and disinfect according to the procedure stipulated in the accompanying Instructions For Use.

This catheter must only be used under medical supervision by personnel who have been adequately trained.

The detailed instructions contained within this IFU is sufficient for guidance of the use of this equipment, no further training from the manufacturer is required. Do not use devices beyond their prescribed lifetime or validated number of uses; injuries may result.

This catheter is to be used in the anatomy and for the age group for which it is intended, as dictated by the device name/description on the label. Intentional use with any other risk of balloon detachment or rupture if balloon filled with water. (Use catheter SR2B if you wish to inflate balloon with water).

TARGET POPULATION

There are no specific intended patient populations. This catheter is for specific patients as at the discretion of the medical professional who have received professional training in gastrointestinal manometry.

The Anorectal-Exulsion Catheter is to be inserted anorectally, for pressure measurements of the rectum. The balloon is to be filled with air to stimulate the rectum.

BALLOON EXPUSION FUNCTION- INSTRUCTIONS FOR USE

- Water-based lubricants may be used to assist in insertion of catheter. Do not use silicone-based lubricants on balloons.
- Insert the catheter through the anus to a depth of 8-10cm.
- Attach a 3-way stopcock to the catheter and fill the balloon with 30-60ml of air (do not fill with water; balloon detachment would occur).
- Close off the balloon port, and remove syringe.
- Have the patient go to the toilet /commode and attempt to expel the catheter.
- Expansion of the catheter within 2-5 minutes indicates normal function.
- Post-procedure manometric catheters are contaminated with body fluid. Immediately clean and disinfect according to the procedure stipulated in the accompanying Instructions For Use.

BIO-FEEDBACK FUNCTION - INSTRUCTIONS FOR USE

For handheld bio-feedback devices:

- With a 3-way stopcock, attach the catheter to the male luer end of the stopcock, a syringe to the female luer end and an extension tube with a male and female luer end to the other female luer end of the stopcock.
- Water-based lubricants may be used to assist in insertion of catheter. Do not use silicone-based lubricants on balloons.
- Insert the catheter through the anus to a depth of 8-10cm.
- Attach the other end of the extension tube to the handheld device and fill the balloon with air accordingly (30-60ml). Close off the syringe port, and remove syringe.
- Do not fill with water; balloon detachment would occur.
- Any squeeze pressure by the patient will be displayed on the handheld device.
- Post-procedure manometric catheters are contaminated with body fluid. Immediately clean and disinfect according to the procedure stipulated in the accompanying Instructions For Use.

CHEMICAL CLEANING / DISINFECTING INSTRUCTIONS

The following reprocessing procedures have been validated by the manufacturer for the user to follow within their own facilities.

- Attach a stopcock to the luer and close it to ensure no liquid enters into the balloon.
- Wash catheter in warm water with enzymatic detergent.
- Rinse catheter with clean fresh water until catheter is visibly clean and disinfected is removed (no residue).
- Wipe external surface.

2. Chemical Disinfection (high-level disinfection)

- Attach the stopcock to the luer and close it to ensure no liquid enters into the balloon.
- Immerse the entire catheter is submerged in disinfectant.
- Soak catheter in disinfectant solution according to solution manufacturer's instructions.

Caution:

- The catheter contains Polyvinyl Chloride (PVC), Polycarbonate, Polysoprene, Silk Suture, Urethane Acrylate Oligomer, N-Dimethylacrylamide and Isobornyl Acrylate.
- Check with disinfectant manufacturer to ensure their disinfectant can be used safely on the materials listed above.
- Reusable PVC Manometric Catheters are compatible for reprocessing with hydrogen peroxide-based disinfectants. Use disinfectants according to the manufacturer's instructions.

DO NOT USE CIDEX PRODUCTS.

- Do not soak the catheter in alcohol.
- Do not steam autoclave the catheter.

3. Final Rinsing and Blow Drying

- Rinse catheter with clean fresh water until catheter is visibly clean and disinfected is removed (no residue).
- Wipe all water from external surface.

STORAGE

- Store catheter in a dry and clean location.

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Limite de température:
15°C - 30°C (60.8°F - 86.6°F)

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Lifetime of product - 5 years or 50 uses whichever occurs first.

Time period for each use - 30 min inside the gastrointestinal system.

It is the responsibility of the user to manually track the age and number of uses of the device (User may contact Mui Scientific for a sampling tracking form).

Manufacturer's Warranty - 6 months from Manufacture Date for any manufacturer's defects.

If user and/or patient should report any serious incident that has occurred in relation to this device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Dispose of the device safely, in accordance with the local policies where this device should be used.

Mui Scientific is not liable for any damage to the catheter, or harm to patients or personnel, caused by improper use of a disinfectant or procedure.

Electronic copy of Instructions For Use is available on Mui Scientific website: www.muiscientific.com

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