EN: REUSABLE PVC MANOMETRIC CATHETERS INSTRUCTIONS FOR USE

INTENDED PURPOSE
Reusable PVC Manometric Catheters are to be used in conjunction with manometric pumps and a computerized data processing system, to measure muscular contractions along the gastrointestinal system. This manometric catheter is designed to be reusable, and to be used only under the supervision of a physician who has received adequate training in gastrointestinal manometry.

USE OF THE REUSABLE PVC
Use of the Reusable PVC Manometric Catheters is indicated when measurement of gastro-intestinal tract pressure are judged to be useful for determining management of patients with proven or suspected gastro-intestinal motor disorders.

CONTRAINDICATIONS
Where there is a significantly increased risk of pulmonary aspiration associated with intubation which cannot be adequately reduced by practical measures such as a longer than usual period of fasting prior to intubation. The presence of any stenosis, structural deformity or disease that results in an unacceptably increased risk of aspiration, perforation or any other adverse consequence arising from use of the device.

WARNINGS
Read through entire IFU to reduce any possible risks from misuse. This device is supplied non-sterile.
Inspect catheter and packaging for any evidence of mechanical damage or imperfections. Do not use if damaged.
The catheter contains low amounts of phthalates. Use with caution for pregnant women.
Reusable PVC Manometric Catheters must not be altered or modified in any way. Mui Scientific is not liable for personal injury and damage to property if original Mui Scientific parts are not being used.
Self-replacement of the catheter balloon is performed at your own risk.
Proper aseptic technique and universal barrier precautions (UBP) must apply.
Post-procedure manometric catheters are contaminated with body fluid. Immediately clean and disinfect according to the procedure stipulated in the accompanying Instructions For Use.

This catheter must only be used under medical supervision by personnel who have been adequately trained in safe gastro-intestinal intubation, and in the performance of perfused manometric measurements.

The detailed instructions contained within this IFU is sufficient for guidance of the catheter and its equipment, not the manufacturer's instructions.
Do not use devices beyond their prescribed lifetime or validated number of uses; injuries may result.
Catheter is to be used in the anatomy and for the age group for which it is intended. Do not use for the device as described on the label.

TARGET POPULATION
There are no specific intended patient populations.
The use for specific patients is at the discretion of the medical personnel who have received professional training in gastrointestinal manometry.

CLINICAL BENEFIT
Manometric Catheters are an integral part of a motility system, for the measurement and assessment of the motor function of the gastrointestinal tract in both, adults and children. The clinical benefit of the Reusable PVC Manometric Catheters, and this overall technique, is to enable the characterization of gastrointestinal motor patterns in the healthy and diseased patients, thus allowing the detection of motility disorders and the assessment of the therapy for patients with problems along their gastrointestinal tract.

INSTRUCTIONS FOR USE
Catheters are to be connected to pressure transducers, which are to be connected to a manometric pump and a computerized data processing system. Together as a system, these devices are used to measure muscular contractions along the gastrointestinal system. The type of pressure transducers and manometric pump used are dictated by the computerized data processing system and the motility study.

• Reusable PVC Manometric Catheters are compatible with any pressure transducer that has a standard luer lock connector. Labeled female luer connector ends of the catheter are to be connected to the corresponding pressure transducer's male luer. A male/female luer connection adaptor may be used, if required.

• Ensure the region of the gastrointestinal tract to be assessed corresponds with the variant type of the catheter device (esophageal, gastric, and small bowel catheters) to be intubated through the nasal or oral cavity, anorectal and colon. Catheters are to be intubated through the rectum.

• Water-based lubricants may be used as insert in insertion of catheter. Do not use silicone-based lubricants on balloons.

• For anorectal catheters with a balloon, the balloon can be filled with air only, to avoid the rectal sensations. The balloon may also be expelled by the patient, to assess the neuromuscular function of the rectum.

• Post-procedure catheters are contaminated with body fluid; only re-use catheters that have been cleaned or reprocessed according to instructions below.

• Do not use silicone-based lubricants on catheters.
• Recommended perfusion flow rate for catheters with 1-12 channels is 0.6mL/min. Do not exceed flow rates greater than 6mL/min.

• The recommended perfusion flow rate for catheters with 13+ channels is 1.5mL/min. Do not exceed flow rates greater than 1.5mL/min.

• Polysoprene (compliant) balloons – inflate up to rated volume; do not exceed double the rated volume; risk of balloon bursting.
• Barostat (non-compliant) balloons – inflate up to rated volume; do not exceed rated volume; do not exceed pressures of 120 mmHg; risk of balloon bursting.

CHEMICAL CLEANING / DISINFECTING INSTRUCTIONS
The following reprocessing procedures have been validated by the manufacturer for user to follow within their own facilities.

1. Cleaning
• If catheter has a balloon attached, plug the end of the luer connector to ensure no liquid enters into the balloon.
• Wash catheter in warm water with enzymatic detergent.
• Flush detergent through all channels including unused channels. Remove all foreign matter. Wipe and clean external surface.

• Rinse catheter and flush all channels with clean fresh water to remove detergent. Repeat as needed.
• Blow-dry all channels with air and wipe excess water from external surface.

2. Chemical Disinfection (high-level disinfection)
• Do not use autoclave for reprocessing of unused channels. Ensure the entire catheter is submerged in disinfectant and all channels are filled with disinfectant.
• Soak catheter in disinfectant solution according to solution manufacturer's instructions.

• The catheter contains Polyvinyl Chloride (PVC), Polycarbonate, Stainless Steel, Urethane Acrylate Oligomer, N-Dimethylacrylamide and Isobornyl Acrylate. If it is not listed on the label, it also contains Polysoprene and Silk Suture.
• Check with disinfectant manufacturer to ensure their disinfectant can be used safely on the materials listed above.

• Reusable PVC Manometric Catheters are compatible for reprocessing with hydrogen peroxide-based disinfectants. Use disinfectants according to the manufacturer's instructions.

DO NOT USE CIDEX PRODUCTS TO DISINFECT ANY CATHETER THAT HAS A BALLOON.
• Do not soak the catheter in alcohol.
• Do not steam autoclave catheters with balloons in open air.

3. Final Rinsing and Blow Drying
• Rinse catheter and flush each channel with sterile water to remove disinfectant. Repeat as needed.
• Blow-dry all channels with air and wipe all water from external surface.

• The following reprocessing procedures have been validated by the manufacturer for the user to follow within their own facilities.

STORAGE
Ensure all channels and external surfaces are completely dry; store catheter dry in open air. Do not store in a sealed container.

Temperature limit:
15°C - 30°C (60.8°F - 86.0°F)

Humidity limit:
10% - 75%

Lifetime of product – 5 years or 50 uses whichever occurs first.
Time period for each use - 30 min inside the gastrointestinal system.
It is the responsibility of the user to manually track the age and number of uses of the device as well as to track the number of uses of the device.

Manufacturer's Warranty – 6 months from Manufacture Date on any manufacturer's defect.
This device and/or patient should report any serious incident that has occurred in relation to this device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
On the date of the device sale, in accordance with the local policies where this device shall be used.
Mui Scientific is not liable for any damage to the catheter, or harm to patients or property caused by use of the device as described on the label.
Electronic copy of Instructions For Use is available on Mui Scientific website: www.muiscientific.com

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