Mui Scientific

Division of
H&A Mui
Enterprises Inc.

145 Traders Blvd. E., Unit #34, Mississauga, ON, Canada, L4Z 3L3
Tel: (905) 890-5525 Toll Free: 1-800-303-6611 Fax: (905) 890-3523
Website: www.muiscientific.com Email: mail@muiscientific.com

EN: SINGLE-USE PVC MANOMETRIC CATHETERS INSTRUCTIONS FOR USE

INTENDED PURPOSE

Single-Use PVC Manometric Catheters are to be used in conjunction with manometric pumps and a computerized data processing system, to measure muscular contractions along the gastrointestinal system. This manometric catheter is designed for single-use, and is to be used only under the supervision of a physician who has received adequate training in gastro-intestinal manometry.

INDICATIONS FOR USE

Use of the Single-Use PVC Manometric Catheters is indicated when measurements of gastrointestinal tract pressures are judged to be useful for determining management of patients with proven or suspected gastrointestinal motor disorder.

CONTRAINDICATIONS

Where there is a significantly increased risk of pulmonary aspiration associated with intubation which cannot be adequately reduced by practical measures such as a longer than usual period of fasting prior to intubation.

The presence of any stenosis, structural deformity or disease that results in an unacceptably increased risk of aspiration, perforation or of any other adverse consequence arising from use of this device.

WARNINGS

Read through entire IFU to reduce any possible risks from misuse.

The device is supplied non-sterile.

Inspect catheter and packaging for any evidence of mechanical damage or imperfections. Do not use if damaged.

The catheter contains low amounts of phthalates. Use with caution for pregnant and pediatric cases.

The catheter must not be altered or modified in any way. Mui Scientific is not liable for personal injury or damage to property if original Mui Scientific parts are not being used.

Apply aseptic technique and universal barrier precautions (UBP) must apply.

The catheter is designed for single use only; do not insert into patient more than once.

Attempts to reprocess, re-sterilize, and/or re-use may lead to failure of the device and/or cross-contamination.

This catheter must only be used under medical supervision by personnel who have been adequately trained in safe gastro-intestinal intubation, and in the performance of perfused manometric measurements.

The detailed instructions contained within this IFU is sufficient for guidance of the usage of this equipment; no further training from the manufacturer is required.

Do not use devices beyond their prescribed lifetime or shelf life; injuries may result.

Catheter is to be used in the anatomy and for the age group for which it is intended as dictated by the device manual/description on the label.

TARGET POPULATION

There are no specific intended patient populations.

The use for specific patients is at the discretion of the medical personnel who have received professional training in gastrointestinal manometry.

CLINICAL BENEFIT

Manometric Catheters are an integral part of motility system, for the measurement and assessment of the motor function of the gastrointestinal tract in both, adults and children. The clinical benefit of the Single-Use PVC Manometric Catheters, and this overall technique, is to enable the characterization of gastrointestinal motor patterns in healthy and diseased patients, thus allowing the detection of motility disorders and the assessment of therapy for patients suffering with problems along their gastrointestinal tract.

INSTRUCTIONS FOR USE

Catheters are to be connected to pressure transducers, which are to be connected to a manometric pump and a computerized data processing system. Together as a system, these devices are used to measure muscular contractions along the gastrointestinal system. The type of pressure transducers and manometric pump used are dictated by the computerized data processing system being used for the motility study. The Single-Use PVC Manometric Catheters are compatible with any pressure transducer that has a standard luer lock connector. Labeled female luer connector ends of the catheter are to be connected to the corresponding pressure transducer's male luer. A male/female luer connector adaptor may be used, if required.

Ensure the region of the gastrointestinal tract to be assessed corresponds with the variant type of the catheter device (esophageal, gastric, and small bowel catheters are to be intubated through the nasal or oral cavity; anorectal and colonic catheters are to be intubated through the rectum).

Do not use silicone-based lubricants on balloons. For anorectal catheters with a balloon, the balloon can be filled with air only, to assess the patient's rectal sensations. The balloon may also be expelled by the patient, to assess the neuromuscular function of the rectum.

Post-procedure catheters are contaminated with body fluid; immediately dispose of device to prevent re-use.

The recommended perfusion flow rate for catheters with 1-12 channels is 0.6mL/min.

The recommended perfusion flow rate for catheters with 13+ channels is 0.15mL/min.

Polyisoprene / Polyolefin (compliant) balloons – inflate up to rated volume; do not exceed double the rated volume; risk of balloon bursting.

STORAGE

Store catheter in a dry and clean location.

Temperature limit: 15°C - 30°C (60.8°F - 86.6°F)

Humidity limit: 10% - 75%

Shelf life: 2 years for catheter with balloon and 3 years for catheter without balloon.

Lifetime: 30 min inside the gastrointestinal system of the human body. Manufacturer's Warranty – 6 months from date of purchase from Mui Scientific for any manufacturer's defects.

The user and/or patient should report any serious incident that has occurred in relation to this device to the manufacturer and the competent authority of the region in which the user and/or patient is established.

Dispose of the device safely, in accordance with the local policies where this device should be used.

Mui Scientific is not liable for any damage to the catheter, or harm to patients or personnel, caused by improper use.

Electronic copy of Instructions For Use is available on Mui Scientific's website: www.muiscientific.com

REF Prefix	Description
SE	Single-Use Esophageal
SR	Single-Use Anorectal
SPE	Single-Use Pediatric Esophageal
SPRB	Single-Use Pediatric Anorectal w/ Balloon
S7	Customized Single-Use

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Authorizes Representative in UK: Advena Ltd., Pure Offices, Plato Close, Tachbrook Park, Warwick, CV34 6WE, UK

Authorizes Representative in EU: Advena Ltd., Tower Business Centre, 2nd Floor, Tower Street, Swatara, BKR 4013, Malta

Authorizes Representative in Switzerland: TAS SAT AG, Chamerstrasse 172, 6300 Zug, Switzerland

Authorizes Representative in USA: Advena Ltd., Tower Business Centre, 2nd Floor, Tower Street, Swatara, BKR 4013, Malta

Authorizes Representative in Canada: Advena Ltd., Tower Business Centre, 2nd Floor, Tower Street, Swatara, BKR 4013, Malta

Authorizes Representative in Mexico: Advena Ltd., Tower Business Centre, 2nd Floor, Tower Street, Swatara, BKR 4013, Malta

Authorizes Representative in India: Advena Ltd., Tower Business Centre, 2nd Floor, Tower Street, Swatara, BKR 4013, Malta

Authorizes Representative in Australia: Advena Ltd., Tower Business Centre, 2nd Floor, Tower Street, Swatara, BKR 4013, Malta

Authorizes Representative in New Zealand: Advena Ltd., Tower Business Centre, 2nd Floor, Tower Street, Swatara, BKR 4013, Malta

Authorizes Representative in South Africa: Advena Ltd., Tower Business Centre, 2nd Floor, Tower Street, Swatara, BKR 4013, Malta

Authorizes Representative in Brazil: Advena Ltd., Tower Business Centre, 2nd Floor, Tower Street, Swatara, BKR 4013, Malta

Authorizes Representative in Argentina: Advena Ltd., Tower Business Centre, 2nd Floor, Tower Street, Swatara, BKR 4013, Malta

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