

**POLYETHYLENE BAROSTAT BALLOONS**

**SQUARE TYPE BALLOON:**

<u>PART #</u>	<u>DIAMETER (cm)</u>	<u>LENGTH (cm)</u>	<u>VOLUME (ml)</u>
1. CT-B40E	6	6	40
2. CT-B60E	10	9	60
3. CT-B80E	6	10	80
4. CT-B120E	6	11	120
5. CT-B150E	6	13	150
6. CT-B250R	6	15	250
7. CT-B500R	10	15	500
8. CT-B600R	10	15	600
9. CT-B800G	10	23	800
10. CT-B1000G	10	25	1000

**PILLOW TYPE BALLOON:**

<u>PART #</u>	<u>DIAMETER (cm)</u>	<u>LENGTH (cm)</u>	<u>VOLUME (ml)</u>
1. CT-BP500R	15	20	500
2. CT-BP600R	15	22	600
3. CT-BP700R	15	23	700
4. CT-BP800R /G	16	25	800
5. CT-BP1000G	20	25	1000


**E-ESOPHAGEAL**

**R-RECTAL**

**G-GASTRIC**


**POLYETHYLENE BAROSTAT BALLOON**

Square Type



40-1000 ml

Pillow Type



40-1000 ml