

POLYOLIFIN BAROSTAT BALLOONS

SQUARE TYPE BALLOON:

	<u>PART #</u>	<u>VOLUME (mL)</u>	<u>DIA. (cm)</u>	<u>LENGTH (cm)</u>
1.	CT-B15E	15	3	5
2.	CT-B30E	30	4	7
3.	CT-B40E	40	3	9
4.	CT-B60E	60	5	9
5.	CT-B80E	80	4	10
6.	CT-B100E	100	5	11
7.	CT-B120E	120	6	11
8.	CT-B150E	150	6	13
9.	CT-B250R	250	6	14
10.	CT-B500R	500	9	17
11.	CT-B600R	600	10	17.5
12.	CT-B800G	800	9	19
13.	CT-B1000G	1000	9	22


PILLOW TYPE BALLOON:

	<u>PART #</u>	<u>VOLUME (mL)</u>	<u>DIA. (cm)</u>	<u>LENGTH (cm)</u>
1.	CT-BP500R	500	15	20
2.	CT-BP600R	600	15	22
3.	CT-BP700R	700	15	23
4.	CT-BP800R/G	800	16	21.5
5.	CT-BP1000G	1000	17.5	23

“E” - ESOPHAGEAL
 “R” - RECTAL/COLONIC
 “G” - GASTRIC

BAROSTAT BALLOON SHAPES

Square Type



Pillow Type

